

CERTIFICATE OF ANALYSIS

Product

2,4-DICHLOROPHENOXYACETIC ACID

Product number

D0911

Batch number

Date of production

019355

Date of expiry

18-01-2025

Packaging and storage

31-01-2029

store in a well-closed container at room temperature

SPECIFICATION

RESULT

UNITS REMARKS

Description

white to off-white or

complies

Loss on Drying

pale yellow powder ≤ 1.0

0.3

%

Assay

 ≥ 97

97.5

This certificate of analysis is conform the results of the manufacturer.

Conclusion: approved

Released by:

J. van der Eem, MSc, pharmacist

Department of Quality Control Duchefa Biochemie B.V.

Date: 18-06-2025

The buyer is responsible to test each batch to ensure the product is suitable for their specific purpose.

Duchefa Biochemie B.V.

A. Hofmanweg 71 2031 BH Haarlem The Netherlands

Postbus 809 🗵 2003 RV Haarlem The Netherlands