



Certificate of Analysis

1.07299.0000 4-(Methylamino)phenol sulfate Technipur™
Batch K55886999

	Spec. Values		Batch Values	
Assay (cerimetric)	≥ 99	%	100	%
Identity (IR-spectrum)	conforms		conforms	
in ether soluble components	≤ 0.2	%	< 0.2	%
pH-value (5 %; water; 20 °C)	3.5 - 4.5		3.6	
Chloride (Cl)	≤ 0.01	%	< 0.01	%
Heavy metals (as Pb)	≤ 0.001	%	< 0.001	%
Fe (Iron)	≤ 0.005	%	< 0.001	%
p-Aminophenol	≤ 2.5	%	≤ 2.5	%
N , N-Dimethyl-p-phenylen diammonium sulfate	passes test		passes test	
Sulfated ash	≤ 0.1	%	< 0.1	%
Loss on drying	≤ 0.3	%	0.1	%

Date of release (DD.MM.YYYY) 02.04.2024
Minimum shelf life (DD.MM.YYYY) 30.04.2026

Dr. Michael Memmel
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.