
PhytoTechnology Laboratories®

and all the

"Helping To Build A Better Tomorrow Through Plant Science"™

Product Information Sheet

A120 Adenine



Synonym:	6-Aminopurine; Vitamin B ₄
CAS:	73-24-5
Formula:	$C_5H_5N_5$
Molecular Wt:	135.13

Properties

Form: Powder Appearance: White to Off-white Powder Application: Cytokinin Solubility: 1N KOH Typical Working Concentration: Storage Temp: Room Temperature Storage Temp of Stock Solution: Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *Phyto*Technology Laboratories® tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References

Merck **13**, 152 Muoma, J., G. Muluvi and J. Machuka, 2008. *In vitro* regeneration by indirect organogenesis of selected kenyan maize genotypes using shoot apices. *Biotechnology*, 7: 732-738.