



Certificate of Analysis

1.01115.0000 Ammonium acetate EMPLURA®
Batch A1679115

	Spec. Values		Batch Values	
Assay (alkalimetric)	≥ 96.0	%	99.1	%
Identity	passes test		passes test	
pH-value (5 %; water)	6.0 - 7.5		6.6	
Chloride (Cl)	≤ 0.002	%	≤ 0.002	%
Sulfate (SO ₄)	≤ 0.01	%	≤ 0.01	%
Heavy metals (as Pb)	≤ 0.0005	%	≤ 0.0005	%
Fe (Iron)	≤ 0.001	%	≤ 0.001	%
Residue on ignition (as sulfate)	≤ 0.02	%	≤ 0.02	%
Water	≤ 2.5	%	≤ 2.5	%

Date of release (DD.MM.YYYY) 26.11.2020
Minimum shelf life (DD.MM.YYYY) 31.10.2025

Claudia Wiegand
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.