

BIOMARK LABORATORIES, INDIA

CERTIFICATE OF ANALYSIS

PRODUCT : Bolton selective Supplement

PRODUCT CODE NO: BF158

QUALITY CONTROL:

PARAMETERS	SPECIFICATIONS/STANDARD VALUES
COLOUR	Offwhite to cream
APPEARANCE	Powder
CLARITY	NA

QUALITY CONTROL : MICROBIOLOGICAL CULTURAL RESPONSE:

Cultural characteristics observed after after an incubation at 35-37°C for 4-6 hours and then at 41.5°C for 40-48 hours when added to BOLTON BROTH BASE(B1320).

Orgm Desc	Inoculum (CFU)	Growth
Candida albicans (ATCC10231)	>=1000	Inhibited
Campylobacter coli (ATCC 33559)	50-100	Good-Luxuriant
Campylobacter jejuni (ATCC 29428)	50-100	Good-Luxuriant
Escherichia coli (ATCC 25922)	>=1000	Inhibited

This is to certify that this lot passes and it confirms to the above mentioned tests and specifications. The information given here is believed to be correct and accurate, however, both the information and products are offered without warranty for any particulars use, other than that specified in the current Biomark literature. The results reported were obtained at the time of release.

This document has been produced electronically and is valid