

## Certificate of Analysis

Product Name:

Brilliant Cresyl Blue ALD - certified by the Biological Stain Commission

Product Number:

B5388

Batch Number:

SHBS0898

Brand:

SIGMA

CAS Number:

51716-96-2

MDL Number:

MFCD00148901

Formula:


C<sub>17</sub>H<sub>20</sub>ClN<sub>3</sub>O · 0.5 ZnCl<sub>2</sub>

Formula Weight:

385.96 g/mol

Quality Release Date:

09 OCT 2024



Test	Specification	Result
Appearance (Form)	Powder	Powder
Infrared Spectrum	Conforms to Structure	Conforms
Solubility (Color) c = 1mg/ml; Water	Dark Blue to Very Dark Blue	Very Dark Blue
Dye Content	≥ 60 %	74 %
Biological Stain Commission Commission Certified Number: AcV-39	Pass	Pass

Larry Coers, Director  
Quality Control  
Milwaukee, WI US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

