

Certificate of Analysis

Product Name:
Carbol Fuchsin

Product Number: C4165
Batch Number: SHBS4380
Brand: SIGMA
MDL Number: MFCD00143923
Quality Release Date: 11 APR 2025
Recommended Retest Date: APR 2029

Test	Specification	Result
Appearance (Color) Green to Brown	Conforms to Requirements	Very Dark Brown
Appearance (Form)	Powder	Powder
Infrared Spectrum	Conforms to Structure	Conforms
Solubility (Color) Dark Pink to Red c = 1mg/ml; Ethanol	Conforms	Conforms
Recommended Retest Period 4 Years	-----	-----



Larry Coers, Director
Quality Control
Sheboygan Falls, WI US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

