


Product Information Sheet

C242 Cyanocobalamin



Synonym: 5,6-Dimethylbenzimidazolyl Cyanocobalamide; Vitamin B₁₂
CAS: 68-19-9
Formula: C₆₃H₈₈CoN₁₄O₁₄P
Molecular Wt: 1355.36

Properties

Form: Powder
Appearance: Red Crystalline
Application: Vitamin
Solubility: Water
Typical Working Concentration:
Storage Temp: 2 to 6° C
Storage Temp of Stock Solution:
Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck **13**, 10074

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™