

## Certificate of Analysis

Product Name:

D-(-)-Salicin -  $\geq 99\%$  (GC)

Product Number:

S0625

Batch Number:

MKCX9587

Brand:

SIGMA

CAS Number:

138-52-3

MDL Number:

MFCD00006590

Formula:


C<sub>13</sub>H<sub>18</sub>O<sub>7</sub>

Formula Weight:

286.28 g/mol

Quality Release Date:

09 JUN 2025



Test	Specification	Result
Appearance (Color)	White to Off White	White
Appearance (Form)	Powder	Powder
Solubility (Color)	Colorless to Faint Yellow	Colorless
Solubility (Turbidity)	Clear to Slightly Hazy	Clear
50 mg/mL, H <sub>2</sub> O		
Infrared Spectrum	Conforms to Structure	Conforms
Purity (GC)	$\geq 99\%$	100 %



Larry Coers, Director

Quality Control

Milwaukee, WI US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

