


Certificate of Analysis

Product Name : Dextrin from corn Type I, powder
Product Number : D2006-VAR
Batch Number : 0000453642
Source Batch : 0000431377
CAS Number : 9004-53-9
MDL Number : MFCD00081554
Quality Release Date : 05 May 2025



Test	Specification	Result
Appearance (Color)	White to Yellow	Off White
Appearance (Form)	Powder	Powder
Residue on ignition (Ash)	≤ 0.50 %	0.00 %
Loss on Drying	≤ 13.0 %	5.8 %
Reducing Sugars (Glucose)	≤ 5 %	2 %

Jagodige Yasomanee, Supervisor

Quality Assurance

St. Louis, Dekalb

US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchase must determine the suitability of the product for its particular use. See reverse side of website or packing slip for additional terms and conditions of sale

Version Number: 01 Doc: 1222485

Page 1 of 1

