



Certificate of Analysis

1.00983.0000 Ethanol absolute for analysis EMSURE® ACS,ISO,Reag. Ph Eur
Batch I1425383

	Spec. Values		Batch Values	
Purity (GC)	≥ 99.9	%	99.9	%
Identity (IR)	conforms		conforms	
Appearance	conforms		conforms	
Color	≤ 10	Hazen	< 5	Hazen
Solubility in water	conforms		conforms	
Acidity or alkalinity	≤ 30	ppm	≤ 30	ppm
Titration acid	≤ 0.0002	meq/g	0.0001	meq/g
Titration base	≤ 0.0002	meq/g	< 0.0002	meq/g
Density (d 20 °C/20 °C)	0.790 - 0.793		0.791	
UV absorption	conforms		conforms	
Aldehydes (as Acetaldehyd)	≤ 0.001	%	≤ 0.001	%
Fusel oils	conforms		conforms	
Substances reducing potassium permanganate (as O)	≤ 0.0002	%	≤ 0.0002	%
Substances reducing permanganate (ACS)	conforms		conforms	
Carbonyl compounds (as CO)	≤ 0.003	%	≤ 0.003	%
Readily carbonizable substances	conforms		conforms	
Acetone, Isopropyl Alcohol (ACS)	conforms		conforms	
Acetone (GC)	≤ 0.001	%	< 0.001	%
Ethylmethylketone (GC)	≤ 0.02	%	< 0.01	%
Isoamyl alcohol (GC)	≤ 0.05	%	< 0.01	%
2-Propanol (GC)	≤ 0.01	%	< 0.01	%
Higher alcohols (GC)	≤ 0.01	%	< 0.01	%
Volatile impurities (GC) (Acetaldehyde and Acetal)	≤ 10	ppm	< 10	ppm
Volatile impurities (GC) (Benzene)	≤ 2	ppm	< 1	ppm
Volatile impurities (GC) (Methanol)	≤ 100	ppm	< 50	ppm
Volatile impurities (GC) (Total of other impurities)	≤ 300	ppm	< 100	ppm
Volatile impurities (GC) (disregard limit)	≤ 9	ppm	9	ppm
Chloride (Cl)	≤ 0.3	ppm	< 0.1	ppm
Nitrate (NO ₃)	≤ 0.3	ppm	< 0.1	ppm
Phosphate (PO ₄)	≤ 0.3	ppm	< 0.1	ppm
Sulfate (SO ₄)	≤ 0.3	ppm	< 0.1	ppm
Ag (Silver)	≤ 0.000002	%	≤ 0.000002	%
Al (Aluminium)	≤ 0.00005	%	≤ 0.00005	%
As (Arsenic)	≤ 0.000002	%	≤ 0.000002	%
Au (Gold)	≤ 0.000002	%	≤ 0.000002	%
Ba (Barium)	≤ 0.00001	%	≤ 0.00001	%
Be (Beryllium)	≤ 0.000002	%	≤ 0.000002	%
Bi (Bismuth)	≤ 0.000002	%	≤ 0.000002	%
Ca (Calcium)	≤ 0.00005	%	≤ 0.00005	%

Certificate of Analysis

1.00983.0000 Ethanol absolute for analysis EMSURE® ACS,ISO,Reag. Ph Eur
Batch I1425383

Cd (Cadmium)	≤ 0.000005	%	≤ 0.000005	%
Co (Cobalt)	≤ 0.000002	%	≤ 0.000002	%
Cr (Chromium)	≤ 0.000002	%	≤ 0.000002	%
Cu (Copper)	≤ 0.000002	%	≤ 0.000002	%
Fe (Iron)	≤ 0.00001	%	≤ 0.00001	%
Ga (Gallium)	≤ 0.000002	%	≤ 0.000002	%
In (Indium)	≤ 0.000002	%	≤ 0.000002	%
Li (Lithium)	≤ 0.000002	%	≤ 0.000002	%
Mg (Magnesium)	≤ 0.00001	%	≤ 0.00001	%
Mn (Manganese)	≤ 0.000002	%	≤ 0.000002	%
Mo (Molybdenum)	≤ 0.000002	%	≤ 0.000002	%
Ni (Nickel)	≤ 0.000002	%	≤ 0.000002	%
Pb (Lead)	≤ 0.00001	%	≤ 0.00001	%
Pt (Platinum)	≤ 0.000002	%	≤ 0.000002	%
Sb (Antimony)	≤ 0.000002	%	≤ 0.000002	%
Sn (Tin)	≤ 0.00001	%	≤ 0.00001	%
Ti (Titanium)	≤ 0.000002	%	≤ 0.000002	%
Tl (Thallium)	≤ 0.000002	%	≤ 0.000002	%
V (Vanadium)	≤ 0.000002	%	≤ 0.000002	%
Zn (Zinc)	≤ 0.00001	%	≤ 0.00001	%
Zr (Zirconium)	≤ 0.000002	%	≤ 0.000002	%
Evaporation residue	≤ 0.0005	%	< 0.0001	%
Water	≤ 0.1	%	0.1	%

Date of release (DD.MM.YYYY) 02.06.2025
Minimum shelf life (DD.MM.YYYY) 31.05.2030

Jeannette David
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.