




Product Information Sheet



M542 Magnesium Nitrate Hexahydrate

Synonym: Nitromagnesite
CAS: 13446-18-9
Formula: $\text{Mg}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$
Molecular Wt: 256.41

Properties

Form: Powder
Appearance: White Crystalline
Application: Plant Macronutrient
Solubility: Water
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck **13**, 5697

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™