



Certificate of Analysis

1.06012.0000 Methanol anhydrous for analysis (max. 0.003% H₂O)
Batch I1413612

Batch Values		
Purity (GC)	≥ 99.9	%
Identity (IR)	conforms	
Color	≤ 10	Hazen
Acidity	≤ 0.0005	meq/g
Alkalinity	≤ 0.0002	meq/g
Acetone (GC)	≤ 0.05	%
Ethanol (GC)	≤ 0.1	%
Evaporation residue	≤ 10	mg/l
Water	≤ 0.003	%

Date of release (DD.MM.YYYY) 29.04.2025
Minimum shelf life (DD.MM.YYYY) 30.04.2028

Jeannette David
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.