



Certificate of Analysis

1.00400.0000 Molybdic acid about 85% MoO₃ (containing ammonium molybdate)
EMPLURA®

Batch K56587400

	Spec. Values		Batch Values	
Assay (precipitative titration)	≥ 85	%	96	%
Chloride (Cl)	≤ 0.01	%	≤ 0.01	%
Phosphate, Arsenate, Silicate (as PO ₄)	≤ 0.005	%	≤ 0.005	%
Sulfate (SO ₄)	≤ 0.02	%	≤ 0.02	%
Fe (Iron)	≤ 0.005	%	≤ 0.005	%
NH ₄ (Ammonium)	about 6	%	about 6	%

Date of release (DD.MM.YYYY) 20.03.2025

Minimum shelf life (DD.MM.YYYY) 31.03.2030

Kristin Starke

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.