

BIOMARK LABORATORIES, INDIA

CERTIFICATE OF ANALYSIS

PRODUCT : Novobiocin supplement
PRODUCT CODE NO: BF083

QUALITY CONTROL:

PARAMETERS	SPECIFICATIONS/STANDARD VALUES
COLOUR	White
APPEARANCE	Powder
CLARITY	NA

QUALITY CONTROL : MICROBIOLOGICAL CULTURAL RESPONSE:

Cultural characteristics observed after (18) hours at (42)oC when added to Tryptic soya broth(B1440) along with supplement (BF083)

Orgm Desc	Inoculum (CFU)	Growth
Escherichia coli 0157 NCTC 12900	50 -100	Good - Luxuriant
Staphylococcus aureus ATCC 25923	$\geq 10^3$	Inhibited

This is to certify that this lot passes and it confirms to the above mentioned tests and specifications. The information given here is believed to be correct and accurate, however, both the information and products are offered without warranty for any particulars use, other than that specified in the current Biomark literature. The results reported were obtained at the time of release.

This document has been produced electronically and is valid.