



Certificate of Analysis

Product Name : Nicotinamide $\geq 98\%$ (HPLC), powder
Product Number : N3376-VAR
Batch Number : 0000462133
Source Batch : 0000447111
CAS Number : 98-92-0
MDL Number : MFCD00006395
Molecular Formula : $C_6H_6N_2O$
Formula Weight : 122.12 g/mol
Quality Release Date : 30 May 2025



Test	Specification	Result
Appearance (Color)	White	White
Appearance (Form)	Powder	Powder
Solubility (Color)	Colorless	Colorless
Solubility (Turbidity)	Clear	Clear
50 mg/mL, H ₂ O		
¹ H NMR Spectrum	Conforms to Structure	Conforms
Carbon	57.5 - 60.5 %	58.9 %
Nitrogen	22.4 - 23.5 %	22.7 %
Purity (HPLC)	$\geq 98\%$	100 %



Jagodige Yasomanee, Supervisor

Quality Assurance

St. Louis, Dekalb

US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchase must determine the suitability of the product for its particular use. See reverse side of website or packing slip for additional terms and conditions of sale

Version Number: 01 Doc: 1227312

Page 1 of 1

