


Certificate of Analysis

Product Name : Phytase from wheat ≥ 0.01 unit/mg solid
Product Number : P1259-5G
Batch Number : 0000466839
Source Batch : 0000385874
CAS Number : 9001-89-2
MDL Number : MFCD00131910
Storage Temperature : -20°C
Quality Release Date : 08 Nov 2024

Test	Specification	Result
units/mg Solid	≥ 0.01	0.04
Enzymatic Activity Unit Definition: One unit will liberate 1 micromole of inorganic phosphorus from 1.5 MM phytate per minute at pH 5.15 at 55 deg C		



Jagodige Yasomanee, Supervisor

Quality Assurance

St. Louis, Dekalb

US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchase must determine the suitability of the product for its particular use. See reverse side of website or packing slip for additional terms and conditions of sale

Version Number: 01 Doc: 1230091

Page 1 of 1

