




Product Information Sheet



P846
Potassium Phosphate,
Monobasic Anhydrous

Synonym: Monopotassium Phosphate; Potassium Dihydrogen Phosphate
CAS: 7778-77-0
Formula: KH_2PO_4
Molecular Wt: 136.09

Properties

Form: Powder
Appearance: White Crystalline
Application: Plant Macronutrient
Solubility: Water
Typical Working
Concentration:
Storage Temp: Room Temperature
Storage Temp of
Stock Solution:
Other Notes: Plant Tissue Culture Tested

Application Notes

References

Merck **13**, 7744

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™