



Certificate of Analysis

5.89579.0000 Pyridine anhydrous, 99.8 %
Batch I1397779

	Spec. Values		Batch Values	
Purity (GC)	≥ 99.75	%	99.90	%
Identity (IR)	conforms		conforms	
Color	≤ 20	Hazen	≤ 20	Hazen
Evaporation residue	≤ 0.0005	%	0.0003	%
Water*	≤ 0.003	%	0.001	%

*For 100 mL and 250 mL unit sizes, water content may exceed above specification; not to exceed 0.005 %.

Date of release (DD.MM.YYYY) 05.02.2025
Minimum shelf life (DD.MM.YYYY) 28.02.2027

Jeannette David
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.