


Certificate of Analysis

Product Name:

Sodium dihydrogencitrate - 99%

Product Number: 234265
Batch Number: BCCN6233
Brand: ALDRICH
CAS Number: 18996-35-5
Formula: C₆H₇NaO₇
Formula Weight: 214,11 g/mol
Quality Release Date: 05 AUG 2025
Recommended Retest Date: JUL 2031



Test	Specification	Result
Appearance (Color)	White	White
Appearance (Form)	Powder or Crystals or Chunk(s)	Powder
Titration with NaOH	98.5 - 101.5 %	99.7 %
Infrared Spectrum	Conforms to Structure	Conforms

Dr. Reinhold Schwenninger
Quality Assurance
Buchs, Switzerland CH

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

