



Certificate of Analysis

1.06329.0000 Sodium hydrogen carbonate for analysis EMSURE® ACS, Reag. Ph Eur
Batch K56474629

	Spec. Values		Batch Values	
Appearance of solution	passes test		passes test	
Identity	passes test		passes test	
Assay (acidimetric; calculated on dried substance)	99.7 - 100.3	%	100.1	%
Assay (acidimetric)	99.0 - 101.0	%	100.1	%
Insoluble matter	≤ 0.015	%	< 0.015	%
Carbonate (CO ₃)	passes test		passes test	
Chloride (Cl)	≤ 0.002	%	≤ 0.002	%
Phosphate (PO ₄)	≤ 0.001	%	≤ 0.001	%
Sulfate (SO ₄)	≤ 0.0150	%	≤ 0.0150	%
Phosphate and Silicate (as SiO ₂)	≤ 0.005	%	≤ 0.005	%
Sulfur compounds (as SO ₄)	≤ 0.003	%	≤ 0.003	%
Total nitrogen (N)	≤ 0.0005	%	≤ 0.0005	%
Ammonium (NH ₄)	≤ 0.0005	%	≤ 0.0005	%
Heavy metals (ACS)	≤ 0.0005	%	≤ 0.0005	%
As (Arsenic)	≤ 0.0002	%	≤ 0.0002	%
Ca (Calcium)	≤ 0.0100	%	≤ 0.0100	%
Cu (Copper)	≤ 0.0002	%	≤ 0.0002	%
Fe (Iron)	≤ 0.0005	%	≤ 0.0005	%
K (Potassium)	≤ 0.005	%	≤ 0.005	%
Mg (Magnesium)	≤ 0.005	%	≤ 0.005	%
Pb (Lead)	≤ 0.0005	%	≤ 0.0005	%
Iodine reducing matter (as I)	≤ 0.0065	%	≤ 0.0065	%
Loss on drying (Silica gel)	≤ 0.20	%	< 0.01	%

Corresponds to ACS, Reag. Ph Eur

Date of release (DD.MM.YYYY) 28.01.2025
Minimum shelf life (DD.MM.YYYY) 31.01.2030

Kristin Starke
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.