



Certificate of Analysis

1.06549.0000 Sodium nitrite for analysis EMSURE® ACS, Reag. Ph Eur
Batch K56471149

	Spec. Values		Batch Values	
Assay (manganometric)	≥ 99.0	%	99.5	%
Insoluble matter	≤ 0.01	%	≤ 0.01	%
Chloride (Cl)	≤ 0.005	%	≤ 0.005	%
Sulfate (SO ₄)	≤ 0.005	%	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.001	%	≤ 0.001	%
Ca (Calcium)	≤ 0.002	%	≤ 0.002	%
Fe (Iron)	≤ 0.001	%	≤ 0.001	%
K (Potassium)	≤ 0.001	%	≤ 0.001	%

Corresponds to ACS, Reag. Ph Eur

Date of release (DD.MM.YYYY) 31.01.2025
Minimum shelf life (DD.MM.YYYY) 31.01.2030

Kristin Starke
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.