



Certificate of Analysis

1.06649.0000 Sodium sulfate anhydrous for analysis EMSURE® ACS,ISO,Reag. Ph Eur
Batch K141985749

	Spec. Values		Batch Values	
Assay (alkalimetric)	≥ 99.0	%	99.9	%
Assay (alkalimetric, calculated on dried substance)	98.5 - 101.0	%	100.0	%
Identity	passes test		passes test	
Appearance of solution	passes test		passes test	
Insoluble matter	≤ 0.01	%	≤ 0.01	%
Acidity or alkalinity	passes test		passes test	
pH-value (5 %; water; 25 °C)	5.2 - 8.0		5.7	
Chloride (Cl)	≤ 0.001	%	≤ 0.001	%
Phosphate (PO ₄)	≤ 0.001	%	≤ 0.001	%
Total nitrogen (N)	≤ 0.0005	%	≤ 0.0005	%
Heavy metals (ACS)	≤ 0.0005	%	≤ 0.0005	%
Heavy metals (as Pb)	≤ 0.0005	%	≤ 0.0005	%
As (Arsenic)	≤ 0.0001	%	≤ 0.0001	%
Ca (Calcium)	≤ 0.005	%	≤ 0.005	%
Fe (Iron)	≤ 0.0005	%	≤ 0.0005	%
K (Potassium)	≤ 0.002	%	≤ 0.002	%
Mg (Magnesium)	≤ 0.001	%	≤ 0.001	%
Loss on drying (130 °C)	≤ 0.5	%	≤ 0.5	%
Loss on ignition (800 °C)	≤ 0.5	%	≤ 0.5	%

Corresponds to ACS, ISO, Reag. Ph Eur

Date of release (DD.MM.YYYY) 21.05.2025
Minimum shelf life (DD.MM.YYYY) 25.03.2028

Kristin Starke
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.