



Certificate of Analysis

1.01252.0000 Starch soluble GR for analysis ISO
Batch F2332352

	Spec. Values	Batch Values
Appearance (color)	white	white
Appearance (description)	fine powder	fine powder
pH-value (2 %; water)	6.0 - 7.5	7.0
Sensitivity	passes test	passes test
Reducing matter (as maltose)	max 0.7 %	< 0.7 %
Sulfated ash	max 0.4 %	0.2 %
Loss on drying (105 °C, 2 h)	max 10 %	7 %
Suitability as enzyme substrate (for amylases)	passes test	passes test

Iso reagent

Minimum shelf life (DD.MM.YYYY) 30.11.2029

Dr. Timm Schlegelmilch
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.