



Certificate of Analysis

Product Name : Testosterone $\geq 98\%$
Product Number : T1500-1G
Batch Number : 0000452489
Source Batch : SLCR4561
CAS Number : 58-22-0
MDL Number : MFCD00003654
Molecular Formula : $C_{19}H_{28}O_2$
Formula Weight : 288.42 g/mol
Recommended Retest Date : Jun 2026
Quality Release Date : 27 Jun 2023



Test	Specification	Result
Appearance (Color)	White to Off White	White
Appearance (Form)	Powder	Powder
Solubility (Color)	Colorless to Very Faint Yellow	Colorless
Solubility (Turbidity) 100 mg/mL, CHCl ₃	Clear	Clear
¹ H NMR Spectrum	Conforms to Structure	Conforms
¹³ C NMR Spectrum	Conforms to Structure	Conforms
Purity (HPLC)	$\geq 98 \%$	100 %



Jagodige Yasomanee, Supervisor

Quality Assurance

St. Louis, Dekalb

US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchase must determine the suitability of the product for its particular use. See reverse side of website or packing slip for additional terms and conditions of sale

Version Number: 01 Doc: 1223782

Page 1 of 1

