

Certificate of Analysis

Product Name:

Tin(II) chloride - reagent grade, 98%

Product Number: 208256
Batch Number: MKCV6888
Brand: SIGALD
CAS Number: 7772-99-8
MDL Number: MFCD00011241
Formula: Cl₂Sn
Formula Weight: 189.62 g/mol
Quality Release Date: 14 FEB 2024



Test	Specification	Result
Appearance (Color)	White to Beige and Faint Purple to Light Purple	White
Appearance (Form)	Powder or Crystals or Flakes	Crystals
ICP Major Analysis Confirms Tin Component	Confirmed	Confirmed
Purity Titration by EDTA	97.5 - 102.5 %	101.1 %



Larry Coers, Director
Quality Control
Milwaukee, WI US

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

