

T8007
Tryptic Soy Agar

Synonym: Tryptic Soy Media

Properties:

Form: Solid
Appearance: Cream to Tan Powder
Application: Plant Tissue Culture Media
Solubility: Water
Typical Working Concentration: 40.00 g/L
Storage Temp: 2-8 °C
Other Notes: Preparation of concentrated solutions is not recommended as insoluble precipitates May form.

Formula [mg/L]:

Sodium Chloride	500
Tryptone	1500
Soytone	500
Agar	1500

Application Notes:

This is a general-purpose medium for use with the cultivation of variety of microorganisms.

Please Note: It is the sole responsibility of the purchaser to determine the appropriateness of this product.

PhytoTech Labs Inc.

14610 W 106th St. Lenexa, KS 66215

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442
phytotechlab.com © 2019 PhytoTech Labs Inc.