



Product Information Sheet

T7900 Tryptic Soy Broth I

Synonym: TSB

Properties

Form: Powder
Appearance: Cream to Tan Powder
Application: Bacterial Screening, Microbiology & Phytopathology
Solubility: Soluble in Water
Typical Working Concentration: 25.00 g/L
Storage Temp: Room Temperature
Storage Temp of Stock Solution: Preparation of concentrated solutions is not recommended as insoluble precipitates may form.

Formula (g/L)

Formula	(g/L)
Tryptone	15
Peptone from Soy	5
Sodium Chloride	5

Application Notes

This is a general-purpose medium for the cultivation of fastidious microorganisms.

Please Note: It is the sole responsibility of the purchaser to determine the appropriateness of this product.

References

1. Difco Manual, 1984. Pg. 1027-1028.

PhytoTechnology Laboratories®

P.O. Box 12205; Shawnee Mission, KS 66282-2205

Phone: 1-888-749-8682 or 1-913-341-5343; Fax: 1-888-449-8682 or 1-913-341-5442

Web Site: www.phytotechlab.com

© 2014 PhytoTechnology Laboratories®