



# Certificate of Analysis

1.09277.0100 Türk's solution for leukocyte counting  
Batch HX56625977

|   | Spec. Values  |    | Batch Values |    |
|---|---------------|----|--------------|----|
| pH-value  | 2.5 - 3.5     |    | 2.8          |    |
| Absorption maximum $\lambda_{\max}$                   | 582 - 586     | nm | 585          | nm |
| Absorbance A ( $\lambda_{\max}$ ; diluted 1:25; 1 cm) | 0.450 - 0.600 |    | 0.515        |    |
| Suitability for microscopy (Counting chamber)         | passes test   |    | passes test  |    |

Date of release (DD.MM.YYYY) 13.03.2025  
Expiry date (DD.MM.YYYY) 29.02.2028

Dr. Timm Schlegelmilch  
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.