

## Certificate of Analysis

Product Name:

**p-Anisidine - 99%**

Product Number:

**A88255**

Batch Number:

**BCCN2247**

Brand:

**ALDRICH**


Formula:

**C<sub>7</sub>H<sub>9</sub>NO**


Formula Weight:

**123,15 g/mol**

Quality Release Date:

**02 MAY 2025**

Test	Specification	Result
Appearance (Color)	White to Brown and Faint Grey to Grey	Faint Brown
Appearance (Form)	Powder or Solid or Chunk(s)	Chunk(s)
Titration with HClO <sub>4</sub>	98.5 - 101.5 %	99.8 %
Purity (GC)	≥ 98.5 %	99.9 %
Infrared Spectrum	Conforms to Structure	Conforms



Dr. Reinhold Schwenninger

Quality Assurance

Buchs, Switzerland CH

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

